

Laser system software

Particulars. si. \rightarrow Laser status always shown at the bottom of the software control program window.
 \rightarrow Inter-locks System \rightarrow stops laser from being on while door is open.

Single frequency \rightarrow default. Single pulse of set freq.
 \rightarrow Set manually to any predefined frequencies.
 \rightarrow Use set button to lock in freq.
 \rightarrow Trigger connector used is "TRout (laser)".

Pulse width \rightarrow width controlled by the DAC.

- \rightarrow DACs are precalibrated.
- \rightarrow Width given in %.
- \rightarrow Also correlated to light intensity.
- \rightarrow New value, turn off laser, enable DAC, enter & set number.

TCT set-up software

Movement Parameters \rightarrow ~~xx~~ $x \leq 0 \rightarrow$ starting position, $dx \rightarrow$ step size,
 $Nx \rightarrow$ number of steps.

transition
vector.

$$R = M \cdot I + T$$

\rightarrow rotation matrix.

: Use default settings for movement {speed - step}
Use time between steps ≈ 0.15 to ensure stages are really stationary while measuring.

Voltage Source

Can choose whether the device is initialised with (Tant Devices).

Output Bias \rightarrow Set voltage for bias \rightarrow then hit apply bias now.

Mode \rightarrow linear? exponential? Hysteresis?

Inc. Speed \rightarrow Speed at which voltage is ramped up / down.

Inc. Bias \rightarrow defines bias step used during ramp up / ramp down.

Delay after bias \rightarrow sets delay to wait while current is measured.